



Ritter®

MEDICAL

244 Barrier-Free®
Exam Table

Safely Accommodate Your Patients

All patients deserve the same access to and quality of care. The Ritter 244 Barrier-Free® Exam Table is designed to safely and gently lift patients weighing up to 850 pounds, time and time again.



Safe + Efficient Care for **Your Patients**

AN EXAM TABLE THAT WORKS AS HARD AS YOU DO

With an **850-pound patient weight capacity**, the Ritter 244 Barrier-Free® Exam Table is designed to safely lift patients to assist in providing a more thorough exam and better access to the patient.



Barrier-Free **Accessibility**

With an 18-inch barrier-free low height and a wide, flat transfer surface, the Ritter® 244 is designed to allow patients to transfer onto and off of the exam table with minimal assistance from staff. While the 34-inch high seat height allows patients to be raised to a more comfortable working height for the caregiver.



HEIGHT ADJUSTABLE

The table adjusts from 18 to 34 inches in height for better access to the table and the patient.



POWERED POSITIONING

Fully powered height adjustment, back, tilt and angle adjustments can allow for better patient positioning. Safety switches in the base section stop the table from moving if a slight amount of force is applied to the base cover.



WIDE SEAT

A 32-inch wide seat can help improve patient comfort and support while the seamless upholstery design helps simplify cleaning.



HAND CONTROL

With the ergonomic, intuitive hand control, patient positioning can be simple and efficient.

OPTIONAL ACCESSORIES



ADJUSTABLE ARM

The optional adjustable arm system is designed to offer patients a stable platform for accessing the table and a sense of security during positioning. Arms rotate and lock into one of three positions and are easily removable.



FOOT CONTROL

The optional foot control is designed to help make patient positioning quick, effective and always within reach.



CASTERS

Optional casters can support patients weighing up to 650 pounds and allow the table to be moved for routine cleaning. (The table should not be used to move patients.)



SPECIFICATIONS

Patient Weight Capacity:
850 lb (385 kg)
650 lb (294.8 kg) with Casters

Height:
Minimum 18" (45.7 cm)
Maximum 34" (86.4 cm)

Back Support Angle:
Minimum = 0°
Maximum = 65°

Seat Tilt Angle:
Minimum = -15°
Maximum = 15°

Upholstery Dimensions
Seat and Back Area:
32" W (81.3 cm)

Head and Foot Area:
30" W (76.2 cm)

Overall Length:
80" (203.2 cm)

Floor Clearance from Floor
to Bottom of Base:

Leveling Feet Installed:
2.75" (7 cm)

Caster Accessory Installed:
3.5" (8.9 cm)

Paper Roll Holder Capacity:
18", 21" or 24"
(45.7 cm, 53.3 cm or
61.0 cm)

ACCESSORIES

Adjustable arm system
(right and left side arms
included)


Casters

Foot control

3' power cord

Electrical
120 VAC, 2.5 A, 60 Hz

STANDARD UPHOLSTERY COLORS

-  Mist 813
-  Stone 814
-  Restful Path 858
-  Iron Ore 851
-  Lunar Gray 845
-  Obsidian 857
-  Oasis 843
-  Soothing Blue 855
-  Dream 847
-  Mineral 844
-  Healing Waters 848
-  Shaded Garden 853
-  Citrus 849
-  Sandy Retreat 850
-  Curative Copper 846
-  Cranberry 859
-  Dark Linen 856
-  Latte 860
-  Robust Brown 852
-  Deep Earth 854

ULTRAFREE UPHOLSTERY COLORS

-  Stone 865
-  Restful Path 868
-  Obsidian 867
-  Harbor 861
-  Wheatgrass 862
-  Spice 863
-  Cranberry 869
-  Dark Linen 866
-  Latte 870
-  Branch 864



CARB 93129.2 Phase 2 Compliant and TSCA Title V1 Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO

For more information, contact your Midmark dealer or call: 1.800.MIDMARK

Outside the USA call: 1.937.526.3662 or visit our website: midmark.com

© 2021 Midmark Corporation, Miamisburg, Ohio USA
Products subject to improvement changes without notice.
Litho in USA.

007-10359-00 Rev. A1 (3/21)

Ritter®